

# S-Blade

## Laryngoscope blade



Reusable fibre optic Miller blade.

### Specific features

Steam sterilizable  
Connection standard ISO 7376

### Technical data

Size	2
Blade length	125 mm
Total length	155 mm
Weight	82 g
Material	Stainless steel
Light transmission	Integrated optical fibers

Class I Med.Device according to Dir. 93/42/CEE

S-BLADE FIBRE-OPTIC SIZE 2-  
MILLER BLADE

**DG03712E**



<b>CND Classification</b>	Z12021003
<b>Registration number</b>	147199

Rev.0 (12/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)