

## S-Blade

### Laryngoscope blade



Reusable fibre optic Miller blade.

#### Specific features

Steam sterilizable  
Connection standard ISO 7376

#### Technical data

Size	4
Blade length	180 mm
Total length	203 mm
Weight	112 g
Material	Stainless steel
Light transmission	Integrated optical fibers

Class I Med.Device according to Dir. 93/42/CEE

S-BLADE FIBRE-OPTIC SIZE 4-  
MILLER BLADE

**DG03714E**



<b>CND Classification</b>	Z12021003
<b>Registration number</b>	147202

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)