

LSK 2

Paediatric laryngoscope kit



A small container, created to always have at disposition an essential kit in a exceptionally reduced space made possible by meticulous organisation.

A kit, which can also be integrated with other kits which do not have a space specifically for laryngoscopes.

It is made in expansive polyethylene and Lycra, is compact and shockproof and has a zip closure. Inside there are areas for blades and handles of laryngoscopes which are well protected by the use of elastic loops.

TECHNICAL DATA

Material:	EVA
Dimensions:	200 x 70 x h200 mm
Weight empty:	350 g
Handle:	Yes
Elastic inserts:	Yes
Reflective:	Yes

LSK 2 empty

DG00101A

