

## LSK 1

### Laryngoscope kit ten blades



#### LSK Series

To work well is also to work with style. This is why Spencer offers a range of laryngoscope holders created by the unmistakable hand of our designer. Not only are they attractive to look at but above all perfectly efficient and rich with intelligent working solutions. Ideal for professional use, exceptional for their compactness and the protection they offer, they will also give your work an image of quality.

All kits in the LSK Series which contain a conventional light or fibre-optic adult handle are integrated with the innovative and exclusive silicone Grip Up handle holder which guarantees the user an extremely pleasant, firm ergonomic grip. A container that can boast the qualities of lightness, mobility, functionality, resistance and durability. Due to its design and versatility, it can easily be inserted in all the backpacks, bags and cases, answering the call for different uses. It is made in expansive polyethylene and Lycra, is compact and shockproof and has a zip closure. Inside there are areas for the placing of blades and handles for laryngoscopes which are well protected thanks to the use of elastic loops. In order to ensure a correct positioning, there is also a padded flap between the two shells.

LSK 1 contains

- a conventional light adult handle
- a conventional light paediatric handle
- five Macintosh blades
- five Miller blades

Also available in fibre-optic version.

Dimensions: 220 x 70 x h220 mm

[8.66 x 2.75 x h8.66 in]

Weight complete: 1,5 kg [3.30 lbs]

#### TECHNICAL DATA

Kit contains:

DG03113	LSK 1 LARINGOSCOPE KIT 10 TRADITIONAL BLADES
DG00100	LSK 1 EMPTY LARINGOSCOPE CONTAINER 10 BLADES
DG03100	L-HANDLE LARINGOSCOPE HANDLE TRADITIONAL LIGHT ADULT
DG03101	L-HANDLE LARINGOSCOPE HANDLE TRADITIONAL LIGHT PEDIATRIC
DG03103	GRIPPER ERGONOMIC HANDLE GRIP ADULT STANDARD
DG04000	S-BLADE MILLER BLADE SIZE 0
DG04001	S-BLADE MILLER BLADE SIZE 1
DG04002	S-BLADE MILLER BLADE SIZE 2
DG04003	S-BLADE MILLER BLADE SIZE 3
DG04004	S-BLADE MILLER BLADE SIZE 4
DG05000	C-BLADE MACINTOSH BLADE SIZE 0
DG05001	C-BLADE MACINTOSH BLADE SIZE 1
DG05002	C-BLADE MACINTOSH BLADE SIZE 2
DG05003	C-BLADE MACINTOSH BLADE SIZE 3
DG05004	C-BLADE MACINTOSH BLADE SIZE 4

LSK 1 complete DG03113B

