

LSK 2

Professional laryngoscope Kit



Pediatric laryngoscope kit with semi-rigid case designed for professionals in reanimation.

Specific features

- The case offers protection to the content
- Zip closure for an easy access to the content
- Laryngoscope handle with non-slip finish
- Macintosh and Miller conventional light blades
- Requires 2 batteries type AA 1,5V not included

Technical data

Dimensions 310x185x80 mm

Included inside the kit

- Pediatric laryngoscope handle
- Macintosh and Miller blades sz. 0,1,2
- Ergonomic handle
- Materials Steel, Chromed brass, EVA
- Handle diameter 18 mm (without cover)
- Handle length 150 mm (non-slip part 120mm)

Blades connection standard

UNI EN ISO 7376

Med.dev. CE marked according to Dir. 93/42/CEE

LSK 2 COMPLETE - PAEDIATRIC
LARYNGOSCOPE KIT
DG031158B



Rev.1 (05/03/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it