

# LSK 4

## Professional laryngoscope Kit



The image represents one size for type of blades included in the kit

Adult laryngoscope kit with semi-rigid case designed for professionals in reanimation.

### Specific features

- The case offers protection to the content
- Zip closure for an easy access to the content
- Laryngoscope handle with non-slip finish
- Supplied with ergonomic silicone handle
- Macintosh and Miller conventional light blades
- Requires 2 batteries type C 1,5V not included

### Technical data

|                   |               |
|-------------------|---------------|
| Dimensions        | 310x185x80 mm |
| Weight of the kit | 1440 ± 20g    |

#### Included inside the kit

|   |                              |
|---|------------------------------|
| Adult laryngoscope handle                 |                              |
| Macintosh and Miller blades sz. 0,1,2,3,4 |                              |
| Ergonomic handle                          |                              |
| Materials                                 | Steel, Chromed brass, EVA    |
| Handle diameter                           | 30 mm (without cover)        |
| Handle length                             | 160 mm (non-slip part 130mm) |

### Blades connection standard

UNI EN ISO 7376

Med.dev. CE marked according to Dir. 93/42/CEE

LSK 4 COMPLETE - PROFESSIONAL  
LARYNGOSCOPE KIT

**DG03116B**

Registration number



58364

Rev.0 (27/09/2016)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)