

## LSK 3

### Disposable laryngoscope kit



A practical, economical and integral solution to always have at hand a complete set of blades and the handle for laryngoscopes. LSK 3 is a useful and versatile solution which results in a container that is made in a technical material; retentive and super protective. Apart from the clear and intuitive disposition of the pockets, LSK 3 characterizes itself for its compactness and the lightness of the complete set means it can be included in all resuscitation kits. Marked in conformity to UNI EN ISO 7376

#### TECHNICAL DATA

Dimensions	170x60x120 mm
Weight	750g
<b>Equipment</b>	
L-Handle	Disposable fibre-optic adult handle
D-Blade	Disposable fibre-optic Macintosh blades Sz. 2,3 and 4
D-Blade	Disposable fibre-optic Miller blades Sz. 0,1 and 2
Latex free DG03810C	LSK 3 – LARYNGOSCOPE KIT 6 F.OPTIC DISP.BLADES



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>