

LSK 3

Disposable blades kit



Practical solution which integrates a set of disposable fiber optic blades.

Specific features

- The case offers great protection to the content
- Zip closure for an easy access to the content
- Macintosh and Miller fiber optic disposable blades
- The kit is compact and light

Technical data

Dimensions | 170x60x120 mm approx

Included inside the kit

- Adult laryngoscope handle
- D-Blade – Disposable fibre optic Macintosh blades Sz 2,3,4
- D-Blade – Disposable fibre optic Miller blades Sz 0,1,2

Blades connection standard

UNI EN ISO 7376

Med.dev. CE marked according to Dir. 93/42/CEE

LSK 3 COMPLETE - DISPOSABLE
LARYNGOSCOPE KIT
DG03810C



Rev.1 (25/05/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it