

MAGILL

Tube-guiding forceps



Magill forceps useful during intubation procedures of infant patients.

Specific features

- Sterilizable
- Reeded distal part to offer a better grip.
- Satin finish

Technical data

Length	140 mm
Curved part length ¹	60 mm
Weight	30 g
Material	Stainless steel

¹ From distal part to first curving

Class I Med.Device according to Dir. 93/42/CEE

MAGILL NEONATAL - TUBE-
GUIDING FORCEPS
FG00606A



CND Classification	L031312
Registration number	141548

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it