

## DGX 550

### Nasopharyngeal airway KIT



Disposable nasopharyngeal airways designed to be inserted into nasal passageway to secure an open airway.

#### Specific Characteristics

Sterile wendle type  
The kit comes in a hard box

#### Technical data

Material	PVC
Internal diameter of each size	5-5,5-6-6,5-7-7,5 mm

Med.Device CE marked according to Dir. 93/42/CEE

DGX 550 KIT SIX SIZES-  
NASOPHARYNGEAL CANNULA

DG90556K



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)