

## Guedel

### 10 Oropharyngeal cannulas



Guedel cannula designed for resuscitation. Once inserted into the oral cavity, prevents the tongue from falling obstructing the upper airway allowing proper airflow.

#### Specific features

- Color coded for a fast identification of the size
- Rounded terminal part to limit possible tissue damage
- Can be sterilized with sterilizing liquid

#### Technical data

Material	PVC
Length	60 mm
Weight	5 g
Ring colour	Black
Included cannulas	10 pz
Latex Free	

Class I Med.Device according to Dir. 93/42/CEE

GUEDEL – DISP. OROPH.  
AIRWAY 60 mm box. 10Pcs.

**RM10321A**



**CND Classification** R010102

**Registration number** 457916

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)