

## Guedel

### 10 Oropharyngeal cannulas



Guedel cannula designed for resuscitation. Once inserted into the oral cavity, prevents the tongue from falling obstructing the upper airway allowing proper airflow.

#### Specific features

- Color coded for a fast identification of the size
- Rounded terminal part to limit possible tissue damage
- Can be sterilized with sterilizing liquid

#### Technical data

Material	PVC
Length	80 mm
Weight	10 g
Ring colour	Green
Included cannulas	10 pz
Latex Free	

Class I Med.Device according to Dir. 93/42/CEE

GUEDEL – DISP. OROPH.  
AIRWAY 80 MM BOX 10 PCS  
RM10323A



CND Classification	R010102
Registration number	457920

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)