

# Guedel

## 10 Oropharyngeal cannulas



Guedel cannula designed for resuscitation. Once inserted into the oral cavity, prevents the tongue from falling obstructing the upper airway allowing proper airflow.

### Specific features

- Color coded for a fast identification of the size
- Rounded terminal part to limit possible tissue damage
- Can be sterilized with sterilizing liquid

### Technical data

Material	PVC
Length	100 mm
Weight	14 g
Ring colour	Red
Included cannulas	10 pz
Latex Free	

Class I med. Dev according to Dir. 93/42/CEE

GUEDEL - DISP. OROPH.  
AIRWAY 100mm BOX 10PCS  
RM10327A



CND Classification	R010102
Registration number	457923

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)