

LMK

Laryngeal masks kit



Kit of disposable laryngeal masks designed to keep patient's airway open during anesthesia or unconsciousness.

Specific features

Inflatable pillow

Disposable and sterile device

The kit comes in a semirigid case characterized by an internal separator equipped with pockets with elastic bands and strap bands to store accessories. The case has a carrying handle and shoulder strap

Reflective insert on the case to increase visibility

Latex Free

Technical data

Case dimensions	360x100x320 ± 10mm
Weight of the kit	1,2 ± 0,1 kg
Included sizes	1-2-3-4-5
Masks material	PVC
Case material	EVA

Med.Device CE marked according to Dir. 93/42/CEE

LMK - LARYNGEAL MASK KIT

RM10610C

Distributed by Spencer Italia s.r.l.

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it