

# DGX 31

## Intubation Stylet



Intubation stylet used to make the application of endotracheal tubes easier.

### Specific features

- Flexible aluminum core, covered by a plastic coating
- Rounded distal parts to avoid tissue damage
- The device is disposable and is provided sterile

### Technical data

Package dimensions	65x495 ± 5 mm
Weight	9 ± 3 g
Length	34 ± 1 cm
Diameter	ø3,3 mm (10Fr)

Class IIa Med.Device according to Dir. 93/42/CEE

DGX 31 PAEDIATRIC - DEFORMABLE  
STYLET

DG50301A

CND Classification



L031199

Distributed by Spencer Italia S.r.l.

Rev.1 (15/03/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)