

RES-Q-BOARD

Resuscitation board



Device designed to optimize CPR maneuvers.

Specific features

- Rigid and anatomical structure
- Dedicated seat for hyperextension of the head for a more effective approach to CPR
- The integrated handles make it easy to handle and transportable
- The conformation of the surface of support and the rigid construction, they offer an optimal support for the patient

Accessories

RB70010A – RES-Q-BOARD – Wall support

Technical data

| | |
|-----------|-------------------|
| Dimension | 605x415x70 ± 5 mm |
| Weight | 1,49 ± 0,2 kg |
| Material | PE |

MD CE marked according to Dir. 93/42/CEE

RES-Q-BOARD – EMERGENCY
BOARD
RB70000A



| | |
|---------------------|--------|
| CND Classification | V0804 |
| Registration number | 136828 |

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it