


RES-Q-BOARD

Resuscitation board



Device designed to optimize CPR maneuvers.

Caratteristiche specifiche

Rigid and anatomical structure
 Dedicated seat for hyperextension of the head for a more effective approach to CPR
 The integrated handles make it easy to handle and transportable
 The conformation of the surface of support and the rigid construction, they offer an optimal support for the patient

Accessories

RB70010A – RES-Q-BOARD – Wall support

Technical data

Dimension	605x415x70 ± 5 mm
Weight	1,49 ± 0,2 kg
Material	PE

MD CE marked according to Dir. 93/42/CEE

RES-Q-BOARD MILITARY -
 RESUSCITATION BOARD
RB70002A



CND Classification	V0804
Registration number	136829

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it