

EVX 30

Peep valve



PEEP Valve used to increase the volume of gas remaining in the lungs at the end of expiration.

Specific features

Universal connector
360° rotation for a more practical use
Autoclavable

Technical data

Dimensions	145x65x45 mm
Adjustment	60cmH ₂ O
Materials	PC, Silicone
Weight	63 g

Class IIa med. dev. according to Dir. 93/42/CEE

EVX 30 AUTOCLAVABLE –
PEEP VALVE

EV60030C



CND Classification R0280

Registration number 354793

Rev.2 (18/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it