

## EVX 32

### Peep valve



PEEP Valve used to increase the volume of gas remaining in the lungs at the end of expiration.

#### Specific features

Universal connector  
360° rotation for a more practical use  
Disposable

#### Technical data

Dimensions	145x65x45 mm
Adjustment	10cmH <sub>2</sub> O
Materials	PC, Silicone
Weight	63 g

Class IIa med. dev. according to Dir. 93/42/CEE

EVX 32 DISPOSABLE - PEEP  
VALVE

EV60032C

CND Classification R0280

Registration number 354794



Rev.2 (07/11/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)