

# EVX 11

## Non-rebreathing valve



Non-rebreathing valve ideal for use with manual resuscitators.

### Specific features

Connection standard ISO 15/22 mm  
 Silicone membranes  
 Autoclavable

### Technical data

Dimensions	87x50x85 mm
Overpressure valve	60cmH <sub>2</sub> O
Materials	PC, Silicone
Weight	45 g

Class IIa med. dev. according to Dir. 93/42/CEE

EVX 11 ADULT NON-REBR.  
 VALVE AUTOCLAVABLE  
 WITH PRES. LIMIT FOR EXP.  
 EV50011C



CND Classification	R0280
Registration number	202025

Rev.2 (17/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)