

EVX 11

Non-rebreathing valve



Non-rebreathing valve ideal for use with manual resuscitators.

Specific features

Connection standard ISO 15/22 mm
 Silicone membranes
 Autoclavable

Technical data

Dimensions	87x50x85 mm
Overpressure valve	60cmH ₂ O
Materials	PC, Silicone
Weight	45 g

Class IIa med. dev. according to Dir. 93/42/CEE

EVX 11 ADULT NON-REBR.
 VALVE AUTOCLAVABLE
 WITH PRES. LIMIT FOR EXP.
 EV50011C



CND Classification	R0280
Registration number	202025

Rev.1 (06/10/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it