

Rock Pin R

Spinal board



The spinal board Spencer Rock R is a device for recovery, immobilization, and transfer to the stretcher of patient with suspected spinal injuries.

Specific characteristics

The number of handles, their positioning and their shape, allow a firm grip even in difficult conditions

The tapered shape allows to save space and to place the device inside Spencer basket stretchers

Dedicated holes and pins to fix belts and Spencer head immobilizers

Smooth and filleted surfaces to allow easy cleaning

Robust and lightweight structure to ensure impact resistance

This model adopts new manufacturing solutions for increased rigidity

Optionals

ST02105A	10g vertical fastener for spinal boards
ST02108A	10g horizontal fastener for spinal boards
ST02015A	PIN STRAPS- SYSTEM OF BELTS WITH HOOKS
ST02035A	Rock straps - straps system for spinal boards
SH00200A	Spencer contour – head immobilizer
ST02106A	Mark Up – system for perm.personalisation

Technical data

Dimensions	1840x445x50 ± 10 mm
Materials	PE, PU
N° of handles	16
N° of pins	14
Loading capacity	200 kg
Weight	8 kg

Class I Med.Device according to Dir. 93/42/CEE

ROCK PIN R - MOD.Z001
REINFORCED YELLOW SPINAL
BOARD

ST02009B

CND Classification



V0804

Rev.0 (13/07/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it