


B-BAK PIN MAX

Spinal board with Pin



The spinal board Spencer B-Bak Pin MAX is a device for recovery, immobilization and transfer to a main stretcher of patients with suspected spinal injuries. This particular version is characterized by increased loading capacity

Specific features

- The number of handles, their placing and their shape offer an easy grip even in difficult conditions
- The tapered shape allows to place the device inside Spencer basket stretchers
- Smooth and rounded surfaces ensure an effective sanitation
- Robust and light structure to ensure high resistance
- Slides on the lower side to minimize sliding resistance during the application
- The low thickness makes the application and storage easier
- Equipped with Pins to allow easy fixing of belts with snap hook

Technical data

Dimensions	1840x405x45 mm
Materials	PE
N° Handles	14
N° PIN	8
Loading capacity	454 kg
Weight	6 ± 0,3 kg

Optional accessories

ST02105A	Fix board - 10g vert.fastener for sp. boards
ST02108A	Fix board - 10g horiz. Fastener for sp.boards
ST02035A	Rock straps – Universal belt system
ST00494B	STX 494 - 3 belts with snap hook
ST02106A	Mark-Up – Perm. Personalization system
SH00200A	Spencer contour – universal head immobilizer

Class I Med.Device according to Dir. 93/42/CEE

B-BAK PIN MAX YELLOW SPINE BOARD

ST02067B

CND Classification



V0804

Rev.2 (25/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it