


## B-BAK PIN MAX

Spinal board with Pin



The spinal board Spencer B-Bak Pin Max is a device for recovery, immobilization and transfer to a main stretcher of patients with suspected spinal injuries. This version is characterized by increased loading capacity if compared with the standard model.

### Specific features

The number of handles, their placing and their shape offer an easy grip even in difficult conditions

The tapered shape allows place the device inside Spencer basket stretchers

Smooth and rounded surfaces ensure an effective sanitation

Robust and light structure to ensure high resistance

Slides on the lower side to minimize sliding resistance during the application

The low thickness makes the application and storage easier

Equipped with Pins to allow easy fixing of belts with snap hook

### Technical data

Dimensions	1840x405x45 mm
Materials	PE
N° Handles	14
N° PIN	8
Loading capacity	454 kg
Weight	6 ± 0,3 kg

### Optional accessories

ST02105A	Fix board - 10g vert.fastener for sp. boards
ST02108A	Fix board - 10g horiz. Fastener for sp.boards
ST02035A	Rock straps – Universal belt system
ST00494B	STX 494 - 3 belts with snap hook
ST02106A	Mark-Up – Perm. Personalization system
SH00200A	Spencer contour – universal head immobilizer

Class I Med.Device according to Dir. 93/42/CEE

B-BAK PIN MAX MILITARY – SPINAL BOARD WITH PIN

**ST02068B**

**CND Classification**



V0804

Rev.2 (25/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[spencer.it](http://spencer.it)