

## B-BAK PIN MAX

### Spinal board with pins



B-bak Pin Max spine boards are devices designed to lift and immobilize patients with suspected spinal injuries.

#### Specific features

- The number of handles, their placing and their shape offer an easy grip even in difficult conditions
- The tapered shape allows place the device inside Spencer basket stretchers
- Smooth and rounded surfaces ensure an effective sanitation
- Robust and light structure to ensure high resistance
- Slides on the lower side to minimize sliding resistance during the application
- The low thickness makes the application and storage easier
- Equipped with pins to allow easy fixing of belts with snap hook

#### Technical data

Dimensions	1840 x 405 x 45 mm
Materials	PE
Number of handles	14
Number of pins	8
Load capacity	454 kg
Weight	6 kg

Class I MD compliant with UE Reg. 2017/745

B-BAK PIN MAX MILITARY - SPINE  
BOARD WITH PIN

ST02068B

CND Classification



V0804

#### Accessories

ST02105A	Fix board – 10g vert. fastener for spine boards
ST02108A	Fix board – 10g horiz. fastener for spine boards
ST02035A	Rock Straps – Universal belt system
ST00494B	STX 494 – 3 belts with snap hook
ST02106A	Mark-Up – Permanent personalization system
SH00201A	Spencer Contour – universal head immobilizer

Rev.0 (14/06/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[spencer.it](http://spencer.it)