

Baby Go

Pediatric spine board



Baby Go spine boards are devices designed to lift and immobilize patients with suspected spinal injuries.

Specific features

Pediatric spine board with differential occipital areas to adapt to the anatomical features of pediatric patients

The fixing of the head immobilizer to spinal takes place by means of straps bands and is equipped with chin and forehead belts that improve the immobilization of the head

Technical data

Dimensions	1190 x 320 x 45 mm
Number of handles	10
Weight	3 kg
Load capacity	30 kg
Material	Polyethylene

Class I MD compliant with UE Reg. 2017/745

BABY GO PAEDIATRIC SPINE BOARD
W/OUT HEADIMMOBIL.



ST02141B

CND classification	V0804
Registration identifier	101672

Rev.0 (14/06/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it