

Baby Go

Pediatric spine board



Immobilization device intended for recovery and transfer on the stretcher pediatric patients.

Specific features

Pediatric spine board with differential occipital areas to adapt to the anatomical features of pediatric patients
The fixing of the head immobilizer to spinal takes place by means of straps bands and is equipped with chin and forehead belts that improve the immobilization of the head.

Technical data

Dimension	1195x320x45 ± 5 mm
N° handle	10
Weight	2,2 kg
Material	PE, PVC, Nylon

MD class I CE marked according to Dir. 93/42/CEE

BABY GO PEDIATRIC SPINE BOARD
W/OUT HEAD IMMOBILIZER



ST02141B

CND classification	V0804
Registration identifier	101672

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it