

MARK UP

Permanent personalization system



System made of two opposed shells which are permanently coupled to the spine board keeping in view the customization.

Specific features

Sealed and irreversible coupling
 Full integration with spinal board
 Compatible with spine board Spencer
 ROCK, ROCK PIN, TANGO

Technical data

Custom window dimension	62x62 mm
Material	Nylon, Polycarbonate

MARK UP – SYSTEM FOR PERMANENT PERSONALIZATION
 ST02106A

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it