

EVA

Resuscitation kit



Compact semi-rigid professional backpack equipped with devices most commonly used in resuscitation procedures.

Included Equipment		Specific Features	
BG06010	ADULT RUBBER RESUSC. BAG	Semi-rigid case to provide more protection to the content	
BG06020	PEDIATRIC RUBBER RESUSC. BAG	Backpack with 2 adjustable shoulder straps	
BG06030	INFANT RUBBER RESUSC. BAG	Internal separator with elastic bands	
EV50009	OVERPRESSURE VALVE 60 cmH2O	Mash pocket on the internal side of the front half of the backpack	
EV50008	OVERPRESSURE VALVE 40 cmH2O (2X)	Handle on the upper side for an easy transport	
RM20405	SILICON FACE MASK SZ 5		
RM20403	SILICONE FACE MASK SZ 3		
RM20502	RUBBER FACE MASK SZ 2		
RM20503	RUBBER FACE MASK SZ 1		
RM20504	RUBBER FACE MASK SZ 0		
RM30700	MOUNDER - JAW SPREADER		
RM10300	GUEDEL – OROPHARYNGEAL CANNULA 50mm		
RM10302	GUEDEL – OROPHARYNGEAL CANNULA 70mm		
RM10304	GUEDEL – OROPHARYNGEAL CANNULA 90mm		
		Technical data	
		Dimensions	380x300x200 ± 30 mm
		Materials	EVA, Nylon
		Weight of empty backpack	0,73 ± 0,1 kg
		Medical devices CE marked according to dir. 93/42/CE	
		KIT EVA W/3 ECO B-LIFE BAGS GUEDEL AND MOUTH WEDGE	
		CB00209A	

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it