

FIX BOARD

10G Wall support for spinal board



The image can differ from the actual product



Fastener designed for vertical fixation of Spencer spinal boards.

Specific Features

Extremely robust and lightweight
Easy to use and install
Quick locking and unlocking system

Technical data

Upper part dimensions	570x25x65 ± 5 mm
Lower part dimensions	180x50x55 ± 5 mm
Material	Steel, Nylon, brass

Reference standard

EN 1789

Class I Med.Device according to Dir. 93/42/CEE

FIX BOARD 10 G WALL SUPPORT FOR
SPINAL BOARD

ST02105A

CND Classification V9099

Registration number 99875



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it