

JEMS

Kit 3 one-piece adjustable collars



One-piece adjustable cervical collars to be used in first aid and emergency situations for a single patient. They are suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

Fast application thanks to the strap closure
Adjustable in vertical direction
Holes for better transpiration
Padded border to increase patient comfort
Easy to clean

Accessories

QC19010C - COLLAR BAG 1 - SACCA PER COLLARI MONOV.

Technical data

Included sizes	S, M, L
Bag dimensions	750x250 mm
Materials	PE
Weight of the kit	310 g

Class I Med.Device according to Dir. 93/42/CEE

JEMS KIT 3SIZES WITH BAG-ADJUST.RIGID 1PIECE COLL.

QC50100A



CND Classification Y060312

Registration number 195180

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it