

DALLAS

Semi-rigid one-piece cervical collar




Dallas cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features

- Light and comfortable
- Fast application thanks to the strap closure
- The white color allows an immediate cleanliness visual feedback

Technical data

Dimension A	125 mm
Dimension B	100 mm
Dimension C	485 mm
Dimension D	65 mm
Material	Expanded PE
Weight	13 g



Class I MD compliant with UE Reg. 2017/745

DALLAS SIZE S - SEMI-RIGID ONE-PIECE CERV. COLLAR

QC22302A

CND Classification Y060312

Registration number 92293



Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it