

## DALLAS

### Semi-rigid one-piece cervical collar




Dallas cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

#### Specific features

- Light and comfortable
- Fast application thanks to the strap closure
- The white color allows an immediate cleanliness visual feedback

#### Technical data

Dimension A	130 mm
Dimension B	105 mm
Dimension C	525 mm
Dimension D	68 mm
Material	Expanded PE
Weight	18 g



Class I MD compliant with UE Reg. 2017/745

DALLAS SIZE M - SEMI-RIGID ONE-PIECE CERV. COLLAR

QC22303A

CND Classification Y060312

Registration number 92298



Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)