

DALLAS

One-piece cervical collar




One-piece cervical collar to be used in first aid and emergency situations for a single patient. They are suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

- Light and comfortable
- Fast application thanks to the strap closure
- The white color allows an immediate cleanliness visual feedback

Technical data

Dimension A	137 mm
Dimension B	107 mm
Dimension C	575 mm
Dimension D	73 mm
Material	PE exp
Weight	17 g



Class I Med.Device according to Dir. 93/42/CEE

DALLAS SIZE L - SEMI-RIGID ONE-PIECE CERV. COLLAR

QC22304A

CND Classification



Y060312

Registration number

92300

Rev.0 (09/06/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it