

DALLAS

One-piece cervical collar




One-piece cervical collar to be used in first aid and emergency situations for a single patient. They are suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

- Light and comfortable
- Fast application thanks to the strap closure
- The white color allows an immediate cleanliness visual feedback

Technical data

Dimension A	160 mm
Dimension B	130 mm
Dimension C	640 mm
Dimension D	88 mm
Material	PE exp
Weight	33 g



Class I Med.Device according to Dir. 93/42/CEE

DALLAS SIZE XL - SEMI-RIGID ONE-PIECE CERV. COLLAR

QC22305A

CND Classification



Y060312

Registration number

92303

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it