

## DALLAS

### Kit 5 one-piece cervical collars



One-piece cervical collars to be used in first aid and emergency situations for a single patient. They are suitable for the immobilization of the cervical spine in trauma patients.

#### Specific Characteristics

Light and comfortable  
Fast application thanks to the strap closure  
The white color allows an immediate cleanliness visual feedback

#### Accessories

QC19010C - COLLAR BAG 1 - SACCA PER COLLARI MONOV.

#### Technical data

|                   |                 |
|-------------------|-----------------|
| Included sizes    | XS, S, M, L, XL |
| Materials         | PE exp          |
| Weight of the kit | 200 g           |

#### Class I Med.Device according to Dir. 93/42/CEE

DALLAS KIT FIVE SIZES-SEMIRIGID  
1PIECE CERV.COLL.

QC22370A

CND Classification Y060312

Registration number 195188



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)