


BLITZ

The rescuer's backpack



Backpack characterized by reduced thickness and ideal to accommodate devices used by rescuers in first aid. This version is equipped with medical devices for adult patients.

Included Equipment		Technical data	
BG02007A	PRO B-LIFE - PVC AD. RES. W/OVERP.VALVE R-VALVE	Dimensions	480x300x120 ± 20 mm
DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS	Volume	12 ± 1L
DR10365K	NON STERILE LATEX GLOVE SZ. M	Material	Nylon, TPU
GM02101K	GAUZE COMPRESSES STERILE 10X10CM 2 SHEETS	Weight empty	0,8 ± 0,1 kg
GM02215K	REEL OF PLASTER 1,25 cm x 5 m	Weight with equipment	2,28 ± 0,1 kg
IF00100C	T.I.V. SPENCER AUTOMATIC TOURNIQUETTE		
IF00426A	SPENCER ICE PACK PE 3 PCS		
RM10321A	GUEDEL SIZE 60 MM, OROPHARYNGEAL CANN.		
RM10322A	GUEDEL SIZE 70 MM, OROPHARYNGEAL CANN.		
RM10323A	GUEDEL SIZE 80 MM, OROPHARYNGEAL CANN.		
RM10892A	CUFFED MASK SIZE. 3 PVC WITHOUT VALVE		
RM10894A	CUFFED MASK SIZE.5 PVC WITHOUT VALVE		
RM30700A	MOUNDER - JAW SPREADER		
ST00630A	HT ADULT - ISOTHERMAL SHEET		

Specific Characteristics	Medical devices CE marked according to dir. 93/42/CE
Dual transport possibility thanks to the padded shoulder straps and carrying handle	BLITZ - ZAINETTO PRIMO SOCCORSO COMPLETO ADULTO
Adjustable shoulder straps	CB04080K
One capacious front pocket with zip closure	
In the main body of the backpack are present:	
- a mesh pocket on the frontal flap	
- 3 adjustable strap bands and elastic bands to store accessories	
The backpack can be totally opened, to allow an easy access to the content	

Rev.1 (14/11/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it