


BLITZ

The rescuer's backpack



Backpack characterized by reduced thickness and ideal to accommodate devices used by rescuers in first aid. This version is equipped with medical devices for adult patients.

Included Equipment

BG02007A	PRO B-LIFE - PVC AD. RES. W/OVERP.VALVE R-VALVE
DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS
DR10365K	NON STERILE LATEX GLOVE SZ. M
GM02101K	GAUZE COMPRESSES STERILE 10X10CM 2 SHEETS
GM02215K	REEL OF PLASTER 1,25 cm x 5 m
IF00100C	T.I.V. SPENCER AUTOMATIC TOURNIQUETTE
IF00426A	SPENCER ICE PACK PE 3 PCS
RM10321A	GUEDEL SIZE 60 MM, OROPHARYNGEAL CANN.
RM10322A	GUEDEL SIZE 70 MM, OROPHARYNGEAL CANN.
RM10323A	GUEDEL SIZE 80 MM, OROPHARYNGEAL CANN.
RM10892A	CUFFED MASK SIZE. 3 PVC WITHOUT VALVE
RM10894A	CUFFED MASK SIZE.5 PVC WITHOUT VALVE
RM30700A	MOUNDER - JAW SPREADER
ST00630A	HT ADULT - ISOTHERMAL SHEET

Technical data

Dimensions	480x300x120 ± 20 mm
Volume	12 ± 1L
Material	Nylon, TPU
Weight empty	0,8 ± 0,1 kg
Weight with equipment	2,28 ± 0,1 kg

Specific Characteristics

Dual transport possibility thanks to the padded shoulder straps and carrying handle

Adjustable shoulder straps

One capacious front pocket with zip closure

In the main body of the backpack are present:

- a mesh pocket on the frontal flap
- 3 adjustable strap bands and elastic bands to store accessories

The backpack can be totally opened, to allow an easy access to the content

Medical devices CE marked according to dir. 93/42/CE

BLITZ - ZAINETTO PRIMO SOCCORSO COMPLETO ADULTO
CB04080K

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it