


## BLITZ

### The rescuer's backpack



Backpack characterized by reduced thickness and ideal to accommodate devices used by rescuers in first aid.

#### Specific Characteristics

Dual transport possibility thanks to the padded shoulder straps and carrying handle  
 Adjustable shoulder straps  
 One capacious front pocket with zip closure  
 In the main body of the backpack are present:  
 - a mesh pocket on the frontal flap  
 - 3 adjustable strap bands and elastic bands to store accessories  
 The backpack can be totally opened, to allow an easy access to the content

#### Technical data

Dimensions	480x300x120 ± 20 mm
Volume	12 ± 1L
Material	Nylon, TPU
Weight empty	0,8 ± 0,1 kg

BLITZ EMPTY - THE RESCUER'S BACKPACK  
 CB04081K

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)