

# CERVER

## One-piece cervical collar



One-piece cervical collar to be used in first aid and emergency situations. They are suitable for the immobilization of the cervical spine in trauma patients.

### Specific Characteristics


- Light and comfortable
- Fast application thanks to the strap closure
- Access for tracheotomy
- Chin support obtained from a flexible part of the collar
- Easy to clean

### Accessories

- QC19010C - COLLAR BAG 1 – Bag for one piece collars

### Technical data

|                 |                 |
|-----------------|-----------------|
| Size            | XS              |
| Materials       | Vinyl coated PE |
| Height <b>A</b> | 78 mm           |
| Length <b>B</b> | 545 mm          |
| Height <b>C</b> | 135 mm          |
| Weight          | 122 g           |



Class I Med.Device according to Dir. 93/42/CEE

CERVER COLLAR SZ.XS  
QC10100K



|                     |         |
|---------------------|---------|
| CND Classification  | Y060312 |
| Registration number | 100563  |

Rev.0 (05/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)