

# CERVER

## One-piece cervical collar



One-piece cervical collar to be used in first aid and emergency situations. They are suitable for the immobilization of the cervical spine in trauma patients.

### Specific Characteristics


- Light and comfortable
- Fast application thanks to the strap closure
- Access for tracheotomy
- Chin support obtained from a flexible part of the collar
- Easy to clean

### Accessories

- QC19010C - COLLAR BAG 1 – Bag for one piece collars

### Technical data

Size	M
Materials	Vinyl coated PE
Height <b>A</b>	105 mm
Length <b>B</b>	635 mm
Height <b>C</b>	160 mm
Weight	170 g



Class I Med.Device according to Dir. 93/42/CEE

CERVER COLLAR SZ.M

QC10102K



CND Classification Y060312

Registration number 100566

Rev.0 (04/03/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)