

CERVER

One-piece cervical collar



One-piece cervical collar to be used in first aid and emergency situations. They are suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics


- Light and comfortable
- Fast application thanks to the strap closure
- Access for tracheotomy
- Chin support obtained from a flexible part of the collar
- Easy to clean

Accessories

QC19010C - COLLAR BAG 1 – Bag for one piece collars

Technical data

Size	L
Materials	Vinyl coated PE
Height A	102 mm
Length B	625 mm
Height C	150 mm
Weight	172 g



Class I Med.Device according to Dir. 93/42/CEE

CERVER - COLLARE CERVICALE
MONOVALVA MIS.L
QC10103K



CND Classification	Y060312
Registration number	100570

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it