

CERVER

One-piece cervical collar



One-piece cervical collars to be used in first aid and emergency situations. They are suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

Light and comfortable
 Fast application thanks to the strap closure
 Access for tracheotomy
 Chin support obtained from a flexible part of the collar
 Easy to clean

Technical data

Included sizes	XS, S, M, L
Bag dimensions	750x250 mm
Materials	Vinyl coated PE
Weight of the kit	700 ± 20 g

Standard equipment

QC19010C - COLLAR BAG 1 – Bag for one piece collars

Class I Med.Device according to Dir. 93/42/CEE

CERVER KIT FOUR SIZES WITH BAG
 QC10104K



CND Classification	Y060312
Registration number	100576

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it