

# CERVER

## One-piece cervical collar



One-piece cervical collars to be used in first aid and emergency situations. They are suitable for the immobilization of the cervical spine in trauma patients.

### Specific Characteristics

- Light and comfortable
- Fast application thanks to the strap closure
- Access for tracheotomy
- Chin support obtained from a flexible part of the collar
- Easy to clean

### Standard equipment

QC19010C - COLLAR BAG 1 – Bag for one piece collars

### Technical data

Included sizes	XS, S, M, L
Bag dimensions	750x250 mm
Materials	Vinyl coated PE
Weight of the kit	700 ± 20 g

Class I Med.Device according to Dir. 93/42/CEE

CERVER KIT FOUR SIZES WITH BAG  
QC10104K



CND Classification	Y060312
Registration number	100576

Rev.0 (04/03/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)