

SPECTAR

Extrication cervical collar



Bivalve cervical collar designed for use in first aid and emergency situations. It's identified as a solution for the immobilization of the cervical spine in trauma patients.


Specific features

- Anatomical shape of the occipital and mandibular part
- Quick application due to strap closure
- Access for tracheotomy
- Dedicated opening for carotid pulse control
- Comfortable contoured padding
- Color coded for quick identification of the size

Technical data

Size	Stout
Materials	LDPE, PE
Weight	185 g
Dimension A	95 mm
Dimension B	55 mm
Dimension C	125 mm
Width front extended valve	355 mm
Width rear extended valve	355 mm

Note: Dimensions subject to tolerances of ± 10 mm



Accessori opzionali

- JM30028A - Appoggio sternale

Class I Med.Device according to Dir. 93/42/CEE

SPECTAR FOAM - 2 PIECE CERVICAL COLLAR SIZE STOUT

JM30024A

CND Classification Y060312

Registration number 1229067



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it