

# NECK YORK

One-piece cervical collar




One-piece cervical collar to be used in first aid and emergency situations. Suitable for the immobilization of the cervical spine in trauma patients.

### Specific Characteristics

- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding

### Technical data

Size	Pediatric
Materials	PE
Weight	100 g
Dimension A	125 mm
Dimension B	445 mm
Dimension C	95 mm



### Accessories

QC19010C - COLLAR BAG 1 – Bag for one-piece collars

Class I Med.Device according to Dir. 93/42/CEE

NEC YORK SZ.2 - 1 PIECE  
RIGID CERVICAL COLLAR  
PEDIATRIC  
QC19002C



CND Classification	Y060312
Registration number	83084

Rev.0 (05/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)