

## NEC YORK

### Kit 6 one-piece cervical collar



One-piece cervical collar to be used in first aid and emergency situations. Suitable for the immobilization of the cervical spine in trauma patients.

#### Specific features


- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding
- 5 sizes (Infant, pediatric, neckless, 2x short, regular)
- Storage bag

#### Standard Accessories

- QC19010C - COLLAR BAG 1 – Bag for one-piece collars

#### Technical data

Materials	PE , PA
Kit weight (with bag)	760 ± 50gr
Dimension A (±5mm)	125 mm (infant)
	125 mm (Pediatric)
	135 mm (Neckless)
	130 mm (short)
	145 mm (Regular)
Dimension B (±5mm)	450 mm (infant)
	445 mm (Pediatric)
	555 mm (Neckless)
	555 (short)
	555 mm (Regular)
Dimension C (±5mm)	80 mm (infant)
	100 mm (Pediatric)
	95 mm (Neckless)
	80 mm (short)
	115 mm (Regular)



Class I Med.Device according to Dir. 93/42/CEE

NEC YORK KIT 6SIZES WITH BAG-1-PIECE CERV.COLL.

QC19015C

CND Classification



Y060312

Registration number

195215

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)