

## NEC YORK

### One-piece cervical collar




One-piece cervical collar to be used in first aid and emergency situations. Suitable for the immobilization of the cervical spine in trauma patients.

#### Specific Characteristics

- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding

#### Technical data

Size	Neckless
Materials	PE
Weight	120 g
Dimension A	135 mm
Dimension B	555 mm
Dimension C	95 mm



#### Accessories

QC19010C - COLLAR BAG 1 – Bag for one-piece collars

Class I Med.Device according to Dir. 93/42/CEE

NEC YORK SIZE 3 (NECKLESS)  
- 1PIECE CERV. COLLAR



QC19020C

CND Classification

Y060312

Registration number

83086

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)