

NEC YORK

One-piece cervical collar




One-piece cervical collar to be used in first aid and emergency situations. Suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding

Technical data

| | |
|-------------|--------|
| Size | Short |
| Materials | PE |
| Weight | 110 g |
| Dimension A | 130 mm |
| Dimension B | 555 mm |
| Dimension C | 80 mm |



Accessories

QC19010C - COLLAR BAG 1 – Bag for one-piece collars

Class I Med.Device according to Dir. 93/42/CEE

NEC YORK SIZE 4 (SHORT) -
1PIECE CERV. COLLAR



QC19021C

CND Classification Y060312

Registration number 83088

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it