

T- STRAPS

Spider straps system
for spine board



Retention system for the immobilization of patients to be used on spinal boards.

Specific characteristics

- Extremely resistant and robust system
- Strap closures to make rescue operations faster
- Fixed transversal belts
- Adjustable longitudinal belt

Technical data

Number of fixing point	10
Belt width	50 mm
Central belt length	From 560 to 1470 mm
Minimum length "V" belts	270 mm
Maximum length "V" belts	540 mm
Minimum width	540 mm
Maximum width	1150 mm
Material	Nylon
Weight	660 ± 20 g

Note: Dimensions subjected to tolerances of ± 30 mm

Class I Med.Device according to Dir. 93/42/CEE

T-STRAPS - UNIV. RESTR. SYS. FOR
SPINE BOA. IMMOB.

ST02039A

CND classification

V9099

Registration number

99405



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it